


MASTER ABRASIVES

External Cylindrical Grinding Machine

CASE STUDY

#1013

Workpiece:	Flange
Operation:	OD Grinding
Machine:	Micromatic e-Grind 200 CNC
Cycle Time:	28 seconds
Stock Removal:	0.2mm
Hardness:	Dia 1: 55 ~ 60 HRC Dia 2: 18 ~ 25 HRC
CpK:	5.295 (for tolerance of 0.020mm)
Run out W.R.T. Centre:	0.005mm
Surface Roughness (Lead free):	0.46 μ m Ra, R-max 4.33
Dressing Skip:	5 Components



Copyright © 2020 Master Abrasives Ltd

CONTACT US FOR MORE INFORMATION

**MASTER
ABRASIVES**

T:01327-703813

sales@master-abrasives.co.uk

www.master-abrasives.co.uk